

Digital Pregnancy



Digital Pregnancy

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Digital Pregnancy

Preview digital Pregnancy pen

Product Features/Principle?

?Characteristics

1. Over 99% accuracy from the day of your expected period
2. The result could be read at 5 minutes.
3. The result could display as "+" "-" or "Pregnant" "Not pregnant". Unmistakable clear result.
4. Beautiful appearance.

?Principle

Human chorionic gonadotropin(hCG) is a glycoprotein hormone excreted by placenta of the pregnant women which is abound existing in the pregnant women's urine and barely existing in the non-pregnant women's urine. Wondfo one step Preview digital Prengnancy pen is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin in urine by the double antibody sandwich method to diagnose if the woman pregnant. It is a reliable index to assist the clinical pregnancy.