

## Master T



Minitechno Brochure Overview : Laboratory photometer with 32  $\mu$ l flow through cell for End-point, Kinetics, Fixed time and multistandard tests. 200 methods stored. Internal Peltier element pre-adjusted at 25/30/37°C. Equipped with graphic display, printer and 6 positions incubator at 37°C. Touch screen display.

Rating: Not Rated Yet

[Ask a question about this product](#)

### Description

[Minitechno Brochure](#)

#### Overview :

Laboratory photometer with 32  $\mu$ l flow through cell for End-point, Kinetics, Fixed time and multistandard tests. 200 methods stored. Internal Peltier element pre-adjusted at 25/30/37°C.

Equipped with graphic display, printer and 6 positions incubator at 37°C. Touch screen display.

#### Dimensions :

- Width (cm) 42
- Length (cm) 38
- Height (cm) 25
- Weight (Kg) 11
- Power supply (Volts) 110/220