

Rayto RT-7600



Rayto RT-7600 Features: Two counting channels, 23 parameters (including 3 histograms). Large color LCD display with touch screen. Three counting modes: whole blood, pre-diluted blood and peripheral blood. Large memory, multi-format output of test results. Including back flush, soaking and high-volt burning to remove the blockage in the reading chamber. Three doors open design for easy maintenance. External power supply adaptor to decrease electrical noises.

Rating: Not Rated Yet

[Ask a question about this product](#)

Description

Rayto RT-7600

[Rayto RT-7600 Brochure](#)

Features:

- Two counting channels, 23 parameters (including 3 histograms).
- Large color LCD display with touch screen.
- Three counting modes: whole blood, pre-diluted blood and peripheral blood.
- Large memory, multi-format output of test results.
- Including back flush, soaking and high-volt burning to remove the blockage in the reading chamber.
- Three doors open design for easy maintenance.
- External power supply adaptor to decrease electrical noises.