


Compact Epifluorescence MicroImager

Digital Epifluorescence Microscope with integrated Dark Chamber and Touchscreen-PC



Touchscreen control; motorized focusing and xy-table; automatic image acquisition and slide scanning; 10x, 20x or 40x objective; up to 1200x native magnification, 2.3 MPixel CMOS camera, 2 fluorescence channels with long-life LED light-sources (DAPI, FITC, GFP, TRITC, CY3, QDot)

KONTAKT

INFO

FÖRDERUNG



Industriestrasse 8, D99427- Weimar, Germany

Telefon : +49 (0) 3643 85593 12 | Telefon : +49 (0) 3643 85593 79

E-mail: nedal.almohamed@alwex.de | E-mail: p.hofmann@alwex.de

Supported by:



Federal Ministry
for Economic Affairs
and Energy

on the basis of a decision
by the German Bundestag



EUROPA FÜR THÜRINGEN
EUROPÄISCHER FONDS FÜR REGIONALE ENTWICKLUNG