

AUTOstainer



MIRALAB
Lab Equipment

● Product Description

Programmable routine stainer for high flexibility and reliability in histology and cytology. A totally automated system featuring convenience. As soon as the staining procedure is completed, an innovative robot arm transfers the slide racks from the stainer to the washing rack then to the dryer rack and the slide ready for testing. no user intervention is required. Incubation time and agitation flexible rinsing and heated or non-heated fan drying functions, its rinsing well and 4 stain wells allow you to perform a very diverse range of stains

Hardware Specifications

Dimensions (WxDxH)	540 x 430 x 355 mm
Weight	8 Kg
Electrical specifications	AC 220/240 – 50 Hz
Normal operating temperature	15-40 °C
Total slide capacity	20 slides (international sizes)
Reagent capacity	4 reagents
Drying Container	air flow \pm 60°C
Washing Container	200 ml distilled water
Programming Memory	10 programs – 10 steps per program.

